

MediPines Summary of Ainslie Study as Published in CHEST



Title: Validation of a Non-invasive Assessment of Pulmonary Gas Exchange During Exercise in Hypoxia

Published: April 2020 Edition of CHEST

Authors: Connor A. Howe, David B. MacLeod, Liisa Wainman, Samuel J. Oliver, Philip N. Ainslie

Participating Universities: University of British Columbia (Canada), Duke University (United States), Bangor University (United Kingdom)

Background: Pulmonary gas exchange efficiency, determined by the alveolar-to-arterial PO_2 difference ($A-aDO_2$), progressively worsens during exercise at sea-level; this response is further elevated during exercise in hypoxia. Traditionally, pulmonary gas exchange efficiency is assessed through measurements of ventilation and end-tidal gases paired with direct arterial blood gas (ABG) sampling. Since these measures have a number of caveats, particularly invasive blood sampling, the development of new approaches for the non-invasive assessment of pulmonary gas exchange is needed.

Research Question: Is a non-invasive method of assessing pulmonary gas exchange (MediPines AGM100[®]) valid during rest and exercise in acute hypoxia?


Study Design and Methods: Twenty-five healthy participants (10 female) completed a staged maximal exercise test on a cycle ergometer in a hypoxic chamber ($FIO_2=0.11$). Simultaneous ABGs via a radial arterial catheter and non-invasive gas-exchange measurements (AGM100) were obtained in two-minute intervals. Non-invasive gas exchange, termed the O_2 deficit, was calculated from the difference between the end-tidal and the calculated PaO_2 (via pulse oximetry and corrected for the Bohr effect by using the end-tidal PCO_2). Non-invasive O_2 deficit was compared to the traditional alveolar to arterial oxygen difference ($A-aDO_2$) using the traditional Riley analysis.

Results: Under combined conditions of normoxic rest, hypoxic rest and hypoxic exercise, the results revealed strong correlations between the calculated $gPaO_2$ (MediPines AGM100) and directly measured PaO_2 (arterial blood gas).

At hypoxic rest and exercise: strong relationships between MediPines AGM100 ($gPaO_2$) and ABG PaO_2 and O_2 deficit with $A-aDO_2$ remained.

$R^2=0.97$ ($p<0.001$; mean bias =1.70 mmHg; $n = 224$)

Bland-Altman: 0.96 ± 2.75 mmHg



Conclusion Summary: This study found that pulmonary gas exchange efficiency measured using a non-invasive gas exchange monitor provided a valid and reliable measure against directly measured arterial blood gasses during hypoxic exercise. Further, the non-invasive oxygen deficit was strongly correlated with $A-aDO_2$ values obtained from the classic $A-aDO_2$. These results provide promising evidence to support the use of non-invasive gas exchange assessments during hypoxic exercise which may be applicable to both laboratory and clinical patient assessments.