

# MediPines Publication Summary: Ainslie AGM100® Validation Study



**Title:** Validation of a Non-invasive Assessment of Pulmonary Gas Exchange During Exercise in Hypoxia

**Published:** April 2020 Edition of CHEST

**Authors:** Connor A. Howe, David B. MacLeod, Liisa Wainman, Samuel J. Oliver, Philip N. Ainslie

**Participating Universities:** University of British Columbia (Canada), Duke University (United States), Bangor University (United Kingdom)

**Background:** Pulmonary gas exchange efficiency, determined by the alveolar-to-arterial  $PO_2$  difference ( $A-aDO_2$ ), progressively worsens during exercise at sea-level; this response is further elevated during exercise in hypoxia. Traditionally, pulmonary gas exchange efficiency is assessed through measurements of ventilation and end-tidal gases paired with direct arterial blood gas (ABG) sampling. Since these measures have a number of caveats, particularly invasive blood sampling, the development of new approaches for the non-invasive assessment of pulmonary gas exchange is needed.


**Research Question:** Is a non-invasive method of assessing pulmonary gas exchange (MediPines AGM100®) valid during rest and exercise in acute hypoxia?

**Study Design and Methods:** 224 data points were obtained from twenty-five healthy participants who completed a staged maximal exercise test on a cycle ergometer in a hypoxic chamber ( $FIO_2=0.11$ ). Simultaneous ABGs via a radial arterial catheter and non-invasive gas-exchange measurements (AGM100) were obtained in two-minute intervals. Non-invasive gas exchange, termed the  $O_2$  deficit, was calculated from the difference between the end-tidal and the calculated  $PaO_2$  (via pulse oximetry and corrected for the Bohr effect by using the end-tidal  $PCO_2$ ). Non-invasive  $O_2$  deficit was compared to the traditional alveolar to arterial oxygen difference ( $A-aDO_2$ ) using the traditional Riley analysis.

**Results:** Under combined conditions of normoxic rest, hypoxic rest and hypoxic exercise, the results revealed strong correlations between the calculated  $gPaO_2$  (MediPines AGM100) and directly measured  $PaO_2$  (arterial blood gas).

**$R^2=0.97$**   
**(n = 224)**

At hypoxic rest and exercise: strong relationships between MediPines AGM100 ( $gPaO_2$ ) and ABG  $PaO_2$  and  $O_2$  deficit with  $A-aDO_2$  remained.



**Conclusion Summary:** This study found that pulmonary gas exchange efficiency measured using a non-invasive gas exchange monitor provided a valid and reliable measure against directly measured arterial blood gasses at rest and during hypoxic exercise. Further, the non-invasive oxygen deficit was strongly correlated with  $A-aDO_2$  values obtained from the classic invasive approach. These results provide promising evidence to support the use of non-invasive gas exchange assessments which may be applicable to both laboratory and clinical patient assessments.

Clinical Study <https://doi.org/10.1016/j.chest.2020.04.017>